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## Fractionation of Polyethylene by High Temperature High Osmotic Pressure Chromatography

Nobuyuki KAGAWA

A high temperature high osmotic pressure chromatography (HOPC) apparatus has been developed, which can separate high molecular weight compounds into components of narrow molecular weight distribution. A new sample injection unit was devised for the injection of concentrated polymer solutions into the loop. By using the apparatus, the polystyrene and polyethylene samples were successfully fractionated into the respective components of narrow average molecular weight, at room temperature for the former and at elevated temperature of 130 for the latter. The apparatus has proved suitable for the preparative separation of compounds with high molecular weight ; a variety of polymers can be separated satisfactorily into fractions of different average molecular weight, good enough for structural analysis and further processing.

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
Fig. 1 The schematic model of interface between concentrated polymer solution and porous medium

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
Fig. 2 The schematic model in the interior of HOPC column.


Fig. 3 The apparatus for high-temperature high osmotic pressure chromatography.

(a) injection unit, (b) N2 gas cylinder, (c) tube, (d) injection loop, (e) pump, (f) oven, (g) column, (h) washing solvent, (i) block heater, (j) heat controller, (k) injection valve, (l) selection valve, (m) pre-heating loop, (n) waste, (o) fraction collector, (p) filter, (q) degasser, (r) selection valve,

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Table 1 The sample composition of PS - 1 for HOPC

	Mw	(%)
F - 700	6,770,000	5.0
F - 288	2,890,000	25.0
F - 128	1,090,000	25.0
F - 80	710,000	25.0
F - 40	360,000	12.5
F - 20	190,000	7.5

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Table 3 Fractionation conditions for HOPC of PS - 1 solution

Packing material	:	Controlled pore glass CPG - 10
Pore size	:	70 , 100 , 140 , 200nm
Column size	:	length = 300mm , I.D. = 7.8mm
Temperature	:	room temperature
Solvent	:	Tetrahydrofrane ( THF )
Sample concentration	:	200mg/ml
Injection volume	:	1 ml
Flow rate	:	0.2ml/min
Collected volume	:	every 1 ml

Table 4 Fractionation conditions for HOPC of PE - 1 and PE - 2 solutions

Packing material	:	Controlled pore glass CPG - 10
Pore size	:	140nm
Column size	:	length = 300mm , I.D. = 7.8mm
Temperature	:	130
Solvent	:	Xylene
Sample concentration	:	100mg/ml
Injection volume	:	1 ml
Flow rate	:	0.5ml/min
Collected volume	:	every 2.5ml

Table 5 Measurement conditions for SEC of PS - 1

Apparatus	:	CCPD/AS8000/CO8010/UV8010 (TOSOH Co.)
Columns	:	GMH - HR - H × 2 (TOSOH Co.)
Temperature	:	40
UV wavelength	:	254nm
Solvent	:	Tetrahydrofran ( THF )
Sample concentration	:	1 mg/ml
Injection volume	:	50 μl
Flow rate	:	1 ml/min

Table 6 Measurement conditions for SEC of PE - 1 and PE - 2

Apparatus	:	150 - C ALC/GPC (Waters Co.)
Columns	:	GMH - HR - H (S) × 3 (TOSOH Co.)
Temperature	:	140
Solvent	:	1 , 2 , 4 - Trichlorobenzene
Sample concentration	:	0.5mg/ml
Injection volume	:	300 μl
Flow rate	:	0.8ml/min


Fig. 4 Comparison of molecular weight distribution for injected PS - 1 and fractions obtained by HOPC at room temperature, standard polystyrene is shown in parenthesis.

Table 7 Results of fractionation by HOPC used

Sample	Mw	Mn	Mw/Mn
PS - 1 Fr - 1 *	6,360,000	5,300,000	1.2
PE - 1 Fr - 1	2,400,000	1,600,000	1.5
PE - 2 Fr - 1	1,580,000	1,050,000	1.5

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
Fig. 5 Comparison of molecular weight distribution of the first fractions (Fr - 1) of polystyrene ; standard polystyrene is shown in parenthesis.


Fig. 6 Relationship between Mw and Mw/Mn of PS - 1 sample fractionated by HOPC at room temperature.

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
Fig. 7 Comparison of molecular weight distribution of injected PE - 1 and fractions obtained by HOPC at 130°C.

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
Fig. 8 Comparison of molecular weight distribution of injected PE - 2 and the first fraction obtained by HOPC at 130°C.

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