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Development of a New TOYOPEARL Hydrophobic Interaction Chromatographic Resin for Antibody Purification

Emi SAKIMA
Yoshimi HASHIMOTO
Yasutami MITOMA

In view of ever expanding demand for therapeutic antibodies in medical market, we may expect that the chromatographic resins for use in the manufacture and purification of antibodies would also go on increasing in market demand. We have now developed a new grade Toyopearl hydrophobic interaction chromatographic (HIC) resin (TOYOPEARL PPG-600M) good for the use at intermediate stage of antibodies manufacturing process.

The new resin was designed to possess an optimal pore size for antibody molecules and a moderate hydrophobicity that was intermediate between the ether and phenyl-based ligands. Compared with the current Tosoh's HIC resins and competitors' products, this resin exhibited not only higher binding capacity and more efficient recovery but also better selectivity in removing impurities from the antibodies of ascites samples.

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
Fig. 1 Relationship between pore size and effective surface area for binding to proteins.
 a) small size pore
 b) large size pore


Fig. 2 Relationship between -globulins binding capacity and exclusion limit.
 Resin : Polyether bonded TOYOPEARL (All research sample)
 Sample : human -globulins


Fig. 3 Calibration curve of optimized material for antibody binding.

Conditions :
 Column size ; 10.7 mmI.D.*150 mm
 Eluent ; H₂O
 Flow rate ; 0.47 mL/min (linear velocity : 30 cm/h)
 Detection ; RI
 Sample ; Pullulan, ethylene glycol
 Injection ; 40 μL

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
Fig. 4 Separation of mouse ascites fluid on some HIC materials.

Conditions :

Column size ; 7.5 mmI.D.*75 mm

Eluent ; a) 0.1 M Phosphate buffer containing 1.8 M ammonium sulfate (pH 7.0)

b) 0.1 M Phosphate buffer (pH 7.0)

Gradient ; a) to b) 30 min linear

Flow rate ; 1mL/min (linear velocity : 136 cm/h)

Detection ; UV (280 nm)

Sample ; Monoclonal IgG1 (Anti human IgE) in mouse ascites fluid (diluted four times)

Injection ; 100 μ L


Fig. 5 Separation of mouse ascites fluid on some HIC materials.

Conditions :

Column size ; 7.5 mmI.D.*75 mm

Eluent ; a) 0.1 M Phosphate buffer containing 1.8 M ammonium sulfate (pH 7.0)

b) 0.1 M Phosphate buffer (pH 7.0)

Gradient ; a) to b) 30min linear

Flow rate ; 1mL/min (linear velocity : 136 cm/h)

Detection ; UV (280 nm)

Sample ; Monoclonal IgG1 (Anti human IgE) in mouse ascites fluid (diluted four times)

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
Fig. 6 Separation of standard proteins.

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Table 1 Static binding capacity of some HIC resins for human γ -globulins

HIC resin	-Globulins binding capacity (mg/mL-gel)
TOYOPEARL new resin	25
Butyl-650M	34
Phenyl-650M	26
Fractogel EMD Phenyl (S)	17
Phenyl Sepharose 6FF (high sub)	41
Butyl Sepharose 4FF	27
Octyl Sepharose 4FF	17

Adsorption Buffer ; 0.1 M phosphate buffer containing 0.8 M ammonium sulfate (pH7.0)


Fig. 8 Separation of IgG and BSA.
 Conditions :
 Column Size ; 7.5 mmI.D.*75 mm
 Eluent ; a) 0.1M Phosphate buffer containing 1.8 M ammonium sulfate (pH 7.0)
 b) 0.1 M Phosphate buffer (pH 7.0)
 Gradient : a) to b) 60 min linear
 Flow rate ; 1mL/min (linear velocity : 136 cm/h)
 Detection ; UV (280 nm)
 Sample ; BSA 3 mg/nL, monoclonal IgG1 (Anti-LH) 1 mg/mL
 Injection ; 100 μ L

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
Fig. 9 Comparison of breakthrough curves.
 Conditions :
 Column size ; 7.8 mmI.D.*20 cm
 Flow rate ; 2.39 mL/min (linear velocity : 300 cm/h)
 Detection ; UV (280 nm)
 Feed ; monoclonal IgG1 (Anti-LH) 1 mg/mL in 0.1M Phosphate buffer containing 0.8 M ammonium sulfate (pH 7.0)

Table 2 Dynamis binding capacity for monoclonal IgG1 (Anti-LH)

HIC Resin	Dynamic binding capacity (mg/mL-gel, 10%leakage)	Eluted amount (mg/mL-gel,)	Recovery ¹ (%)
TOYOPEARL New Resin	38	39	100
TOYOPEARL Phenyl-650M	30	29	95
Phenyl Sepharose 6FF (high sub)	65	59	89

¹Recovery : (Eluted amount/dynamic binding capacity) × 100


Fig.10 Comparison of elution chromatogram.

Conditions :

Column Size ; 7.8 mmI.D.*20 cm

Flow rate ; 2.39 mL/min (linear velocity : 300 cm/h)

Detection ; UV (280 nm)

Feed ; monoclonal IgG1 (Anti-LH) 1 mg/mL

in 0.1M Phosphate buffer containing 0.8 M ammonium sulfate (pH 7.0)

Wash ; 0.1 M Phosphate buffer containing 1.1 M ammonium sulfate (pH 7.0)

Elute ; 0.1 M Phosphate Buffer (pH7.0)

Feed (until 10%leak) Wash (6 min) Elute (Recovered as shown by arrows)

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