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
Fig. 1 Activation and immobilization reactions

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
Fig. 2 Effect of pH of coupling buffer on (●) antibody coupling yield, (○) antigen binding efficiency and (■) antigen binding capacity. Affinity-purified anti-human albumin antibody (2.5mg) was reacted with 0.25 g of Tressyl-5PW in 2 ml of 1 M potassium phosphate buffer of pH 6.5-8.5 at 4 °C for 16h.


Fig. 3 Structure of IgG

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
Fig. 4 Effect of amount of reacted antibody on (●) antibody coupling yield, (○) antigen binding efficiency and (■) antigen binding capacity. Affinity-purified anti-human albumin antibody (0.5-5mg) was reacted with 0.25g of Tressyl-5PW in 2 ml of 1 M potassium phosphate buffer (pH7.5) at 4 °C for 16h.

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
Fig. 5 Effect of pH of coupling buffer on (●) antibody coupling yield and (○) antigen binding capacity. IgG fraction of anti-human albumin antibody (10.7mg) was reacted with 0.25 g of Tressyl-5PW in 2 ml of 1 M potassium phosphate buffer of pH 6.5-8.5 at 25 °C for 16h.


Fig. 6 Effect of reaction time on (●) antibody coupling yield and (○) antigen binding capacity. IgG fraction of anti-human albumin antibody (10.7mg) was reacted with 0.25 g of Tressyl-5PW in 2 ml of 1 M potassium phosphate buffer (pH 7.5) at 25 °C.

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Table 1 Results of coupling of various antibodies to Tresyl-5PW

Antibody	Amount of antibody reccted [mg]	Coupling yield [%]	Antigen binding capacity [mg/ml gel]
Anti-human albumin antibody	9.0	89	2.8
Anti-human transferrin antibody	7.5	89	2.0
Anti-human IgG antibody	5.5	96	2.5
Anti-human IgM antibody	5.0	90	2.5
Anti-human 1-acidglycoprotein antil	12	83	-
Anti-human plasminogenantibody	5.5	90	-

Anti-human serum protein antibodies (5-12 mg) were reacted with 0.25 of Tresyl-5PW in 1 M potassium phosphate buffer (pH 7.5) at 25 for 8 h. Antigen binding capacity was not determined for anti-human 1-acidglycoprotein antibody and anti-human plasminogen antibody.


Fig. 7 Purification of albumin by high-purance affinity chromatography

Human serum (6 μ l) was applied to a column (6 mm I.D. \times 1 cm) of Tresyl-5PW coupled with anti-human albumin antibody at a flow-rate of 1 ml/min at 25 . The pH of eluent was changed stepwise from 7.0 to 1.6 at 3 min after the sample injection. Bound and unbound peaks were collected and subjected to SDS-PAGE. The results of SDS-PAGE are included; lanes 1 and 5, lane 2, lane 3 and lane 4 represents patterns of molecular weight markers, haman serum, unbound fraction and bound fraction, respectively.


Fig. 8 Purification of transferrin by high-performance affinity chromatography
 Experimental procedure as for albumin in Fig. 7, except that 80 μ l of human serum was applied to a column of Tresyl-5PW coupled with anti-human transferrin antibody. The results of SDS-PAGE are included; lanes 1 and 5, lane 2, lane 3 and lane 4 represents patterns of molecular weight markers, human serum, unbound fraction and bound fraction, respectively.


Fig. 9 Purification of plasminogen by high-performance affinity chromatography
 Experimental procedure as for albumin in Fig. 7, except that 500 μ l of human serum was applied to a column of Tresyl-5PW coupled with anti-human plasminogen antibody and the pH of the eluent was changed at 5 min after the sample injection. The results of SDS-PAGE are included; lanes 1 and 5, lane 2, lane 3 and lane 4 represents patterns of molecular weight markers, human serum, unbound fraction and bound fraction, respectively.

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
Fig.10 Analysis of IgG in human serum by high-performance affinity chromatography
A 1 μ l of human serum diluted 10-fold with 0.1 M phosphate buffer (pH 7.0) was applied to a column (6 mm I.D. \times 1 cm) of Tyesyl-5PW coupled with anti-human IgG antibody at a flow-rate of 1 ml/min at 25 $^{\circ}$ C. Bound IgG was eluted by pulse injection of 1 ml of 0.1 M citric acid (pH 1.6) at 1 min after the sample injection.


Fig.11 Relationship between IgG contents in human serum determined by single radial immunodiffusion (SRID) and by high-performance affinity chromatography (HPAC) Three human serum standards (1 μ l) of different IgG contents were analyzed in the same way as in Fig. 10.

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
Fig.12 SDS-PAGE analysis


Fig.13 Removal of albumin and IgG by high-performance affinity chromatography

Human serum (20 μ l) was applied to a column (6 mm I.D. \times 1 cm) of Tresyl-5PW coupled with anti-human albumin and anti-human IgG antibody at a flow-rate of 0.5 ml/min at 25 $^{\circ}$ C. The pH of eluent was changed stepwise from 7.4 to 1.6 at 8 min after the sample injection. Unbound and bound peaks were collected and subjected to SDS-PAGE.

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