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Determination of Antimony in Chromium Metal by Flameless Atomic Absorption Spectrometry

Takashi HÔJYO
Shinkô NAGUMO
Tokio FUJIKI

A rapid and highly sensitive analytical method has been developed to determine antimony below a few ppm level in electrolytic chromium metal. The method is based on a graphite furnace Zeeman atomic absorption spectrometry and does not require pre-separation of antimony. Preliminary studies using a standard solution containing less than $0.1 \mu\text{g}/\text{ml}$ of antimony revealed a limit of detection to be $0.002 \mu\text{g}/\text{ml}$ under the optimum instrumental operating conditions. When a given sample solution contains more than $5,000 \mu\text{g}/\text{ml}$ of chromium, a modification of drying and ashing conditions is required because of the appearance of bumping and/or too large background absorbance. The detection limit of the present method is 0.4 ppm and the results obtained with several samples are in good agreement with those obtained by a spectrophotometric method (Rhodamine B method).

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
Fig. 1 Effect of heating program

- Drying (30 sec)
- △— Ashing (30 sec)
- Atomizing (7.5 sec)

Table 1 Instrumental operating conditions

	(1) for the standard solution	(2) for the electrolytic chromium metals
Wave length	217.5 nm	
Lamp current	15mA	
Drying	120°C × 30 sec	200°C Ramp mode
Ashing	400°C × 30 sec	650°C × 120 sec
Atomizing	2,540°C × 7.5 sec.	
Argon gas flow		
Sheathing	3 l/min.	
Carrier	0 ml/min.	
Cuvette	Tube type	
Slit width	3	
Expansion	×5	
Response	3	
Sample volume	20 μl	

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
Fig. 3 Effect of chromium

- Chloride } on the condition of Table 1(1)
- △— Nitride } on the condition of Table 1(1)
- Chloride } on the condition of Table 1(2)
- ▲— Nitride } on the condition of Table 1(2)
- ◆— Blank of chloride } on the condition of Table 1(2)
- ▼— Blank of nitride } on the condition of Table 1(2)

Table 2 Analytical results

Sample	Sb added (μg)	Sb found (μg)	Sb (ppm)
Electrolytic chromium metal	0	<0.1	<0.4
	1.0	1.0	
	2.0	2.0	
I	0	<0.1	<0.4
	1.0	1.0	
	2.0	1.8	
II	0	<0.1	<0.4
	1.0	1.0	
	2.0	2.0	
III	0	<0.1	<0.4
	1.0	1.0	
	2.0	2.0	

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