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High-Performance Ion-Exchange Chromatography of Proteins on Non-Porous Ion-Exchanger

Akane MITSUI
Takashi KITAMURA
Yoshio KATO

Porous packing materials have been employed with great success for a long time in liquid chromatography. On the other hand, non-porous packings have been examined and demonstrated to be useful only recently.

Now, we have developed anion and cation exchangers, TSKgel DEAE-NPR and TSKgel SP-NPR, which are prepared by introducing diethylaminoethyl groups or sulfopropyl groups onto non-porous spherical hydrophilic resins of 2.5 μm in diameter. We report herein the characteristics of these new packings and several applications to protein separation.

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
Fig. 1 Titration curve of TSKgel DEAE-NPR


Fig. 2 Titration curve of TSKgel SP-NPR


Fig. 3 Separation of a protein mixture on TSKgel DEAE-NPR

A mixture of conalbumin (1), transferrin (2), ovalbumin (3) and trypsin inhibitor (4) was separated with a 10-min linear gradient of NaCl from 0 to 0.5 M in 20 mM Tris-HCl buffer (pH 8.0) at a flow rate of 1.5 ml/min.

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
Fig. 4 Separation of a protein mixture on TSKgel SP-NPR

A mixture of trypsinogen (1), α -chymotrypsinogen A (2), ribonuclease A (3), cytochrome c (4) and lysozyme (5) was separated with a 10-min linear gradient of sodium sulfate from 0 to 0.5 M in 20 mM acetate buffer (pH 5.0) at a flow rate of 1.5 ml/min.


Fig. 5 Comparison of TSKgel DEAE-NPR and DEAE-5PW

A mixture of ovalbumin (first peak) and trypsin inhibitor (second peak) was separated on DEAE-NPR (left) and DEAE-5PW (right). Elution conditions as in Fig. 3 except the gradient time was 60 min and flow rate was 1.0 ml/min on DEAE-5PW. The size of DEAE-5PW column was 75 \times 7.5 mm I.D.

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
Fig. 6 Comparison of TSKgel SP-NPR and SP-5PW

A mixture of α -chymotrypsinogen A (first peak) and ribonuclease A (second peak) was separated on SP-NPR (left) and SP-5PW (right) with a 10-min (SP-NPR) or 60-min (SP-5PW) linear gradient of NaCl from 0 to 0.5 M in 20 mM 2[N-morpholino]ethanesulfonic acid-NaOH buffer (pH 6.0) at a flow rate of 1.5 ml/min (SP-NPR) or 1.0 ml/min (SP-5PW). The size of SP-5PW column was 75 \times 7.5 mm I.D.


Fig. 7 Dependence of resolution on flow rate in the separation of proteins on TSKgel DEAE-NPR

A mixture of ovalbumin and trypsin inhibitor was separated under the same conditions as in Fig. 3 except the flow rate was varied between 0.25 and 2.0 ml/min. The resolution was calculated from the peak widths and elution volumes of the two proteins.

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
Fig. 8 Dependence of resolution on gradient time in the separation of proteins on TSKgel DEAE-NPR

A mixture of ovalbumin and trypsin inhibitor was separated under the same conditions as in Fig. 3 except the gradient time was varied between 1 and 40 min. Resolutions as in Fig. 7.


Fig. 9 Dependence of resolution on sample loading in the separation of proteins on TSKgel DEAE-NPR

A mixture of equal amount of ovalbumin and trypsin inhibitor was separated under the same conditions as in Fig. 3. The sample loading was varied between 1.25 to 80 μg . Resolutions as in Fig. 7.

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
Fig. 10 Dependence of peak width on sample loading in the separation of proteins on TSKgel DEAE-NPR

Elution conditions as in Fig. 3.

Samples: (○) conalbumin; (●) ovalbumin; (◐) trypsin inhibitor.


Fig. 11 Chromatograms of a crude sample of lipoxidase obtained on TSKgel DEAE-NPR
Sample loadings were 50 μg (left) and 200 μg (right). Elution conditions as in Fig. 3.


Fig. 12 Chromatograms of a crude sample of lipoxidase obtained on TSKgel SP-NPR
A crude sample of lipoxidase was separated on TSKgel SP-NPR with a 10-min linear gradient of sodium sulfate from 0 to 0.5 M in 20 mM acetate buffer (pH 4.5) at a flow rate of 1.5 ml/min. Sample loadings were 50 μg (left) and 200 μg (right).

Table 1 Recovery of proteins from TSKgel DEAE-NPR column

Each protein of 5 μ g was separated on DEAE-NPR under the same conditions as in Fig. 3. The recovery was determined from the areas of eluted peaks.

Protein	Recovery (%)
Thyroglobulin	100
Ferritin	99
γ -Globulin	104
Bovine serum albumin	102
Hemoglobin	91
Ovalbumin	103
β -Lactoglobulin	101
Trypsin inhibitor	98

Table 2 Recovery of proteins from TSKgel SP-NPR column

Each protein of 5 μ g was separated on SP-NPR with a 10-min linear gradient of NaCl from 0 to 0.5 M in 20 mM phosphate buffer (pH 7.0) at a flow rate of 1.5 ml/min. The recovery was determined from the areas of eluted peaks.

Protein	Recovery (%)
Hemoglobin	88
α -Chymotrypsinogen A	95
α -Chymotrypsin	100
Trypsinogen	87
Lysozyme	96
Ribonuclease A	95
Cytochrome c	93


Fig. 13 Separation of protein A from recombinant *E. coli* product (upper) and from *staphylococcus aureus* (lower) by ion-exchange chromatography on TSKgel DEAE-NPR
Conditions as in Fig. 3.

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
Fig. 14 Separation of α -chymotrypsin by ion-exchange chromatography on TSKgel DEAE-NPR

α -Chymotrypsin of 5 μ g was separated with a 20-min linear gradient of NaCl from 0 to 0.5 M in 20 mM 1, 3-diaminopropane-HCl buffer (pH 10.3) at a flow rate of 1.5 ml/min.


Fig. 15 Separation of hexokinase by ion-exchange chromatography on TSKgel DEAE-NPR
A crude sample of hexokinase (25 μ g) was separated under the same conditions as in Fig. 3.

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
Fig. 16 Separation of catalase by ion-exchange chromatography on TSKgel DEAE-NPR
A crude sample of catalase ($50 \mu\text{g}$) was separated under the same conditions as in Fig. 3.


Fig. 17 Purity test of catalase fraction by reversed-phase chromatography on TSKgel Phenyl-5PW RP

Commercial catalase (upper curve) and catalase fraction on TSKgel DEAE-NPR (lower curve) was separated with 2-min linear gradient of acetonitrile from 5 to 20% followed by 48-min linear gradient of acetonitrile from 20 to 80% in 0.05% TFA at a flow rate of 1 ml/min. The size of Phenyl-5PW RP column was 75×4.6 mm I.D.


Fig. 18 Separation of human hemoglobin standard by ion-exchange chromatography on TSKgel SP-NPR

Human hemoglobin standard was separated with a 3-min linear gradient of NaCl from 0.05 to 0.2 M in 20 mM bis-Tris-HCl buffer (pH 6.0) at a flow rate of 1.5 ml/min and detected at 415 nm.

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