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Recent Advance in High-Performance Liquid Chromatography of Proteins, Peptides, and Nucleic Acids on Non-porous Packings

Shigeru NAKATANI
Takashi KITAMURA
Yosuke YAMASAKI
Katsuo KOMIYA
Yoshio KATO

High-performance liquid chromatography on non-porous packings is an efficient tool for the rapid separation of proteins. We have recently demonstrated that the non-porous packings are also useful for the separation of peptides and nucleic acids. In the present paper, new applications of these packings to the separation of proteins, peptides and nucleic acids are reviewed briefly.

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
- (I) Eddy diffusion and mobile-phase mass transfer
 (II) Longitudinal diffusion
 (III) Stagnant mobile-phase mass transfer
 (IV) Stationary-phase mass transfer

Fig. 1. Relationship between HETP and flow rate

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Column size: 4.6 mm I.D. \times 3.5 cm


Fig. 2 Schematic representation of TSKgel NPR Series


Fig. 3 Separation of monoclonal antibody from ascites fluid by high-performance hydroxyapatite chromatography on TSKgel HA-1000

Column : TSKgel HA-1000 (75×7.5 mm I.D.+10×6 mm I.D.)

Sample : anti-chicken 14 K lectin, ascites fluid (1.52 mg in 0.1 ml)

Elution : step gradient of sodium phosphate (pH 5.8) from 60 to 150 mM at 12 min, followed by 300 mM at 24 min and 500 mM at 36 min

Flow rate: 1.0 ml/min

Temperature: 25°C

Detection: UV (280 nm)


Fig. 5 Separation of monoclonal antibody (IgG₁) from ascites fluid by high-performance hydrophobic interaction chromatography on TSKgel Butyl-NPR

Column : TSKgel Butyl-NPR (35×4.6 mm I.D.)

Sample : anti-chicken 14 K lectin, ascites fluid (38 μg in 2.5 μl)

Elution : 10 min linear gradient of ammonium sulfate from 1.5 M to 0 in 0.1 M sodium phosphate (pH 7.0)

Flow rate: 1.0 ml/min

Temperature: 25°C

Detection: UV (280 nm)


Fig. 4 Purity check of monoclonal antibody fraction obtained from TSKgel HA-1000 by high-performance ion-exchange chromatography on TSKgel DEAE-NPR

Column : TSKgel DEAE-NPR (35×4.6 mm I.D.)

Sample : anti-chicken 14 K lectin, A. ascites fluid (76 μg in 5 μl), B. unbound fraction (100 μl), C. bound IgG fraction (100 μl)

Elution : 10 min linear gradient of NaCl from 0 to 0.5 M in 50 mM Tris-HCl (pH 8.5)

Flow rate: 1.5 ml/min

Temperature: 25°C

Detection: UV (280 nm)


Fig. 6 Separation of trypsin by affinity chromatography on TSKgel STI-5PW

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
Fig. 7 Separation of trypsin by hydrophobic interaction chromatography on TSKgel Butyl-NPR

Column : TSKgel Butyl-NPR (35×4.6 mm I.D.)

Sample : bovine pancreas trypsin; (A) original sample (10 μ g), (B) fraction 1, (C) fraction 2

Elution : 10 min linear gradient of ammonium sulfate from 2 to 0 M in 20 mM Tris-HCl buffer, pH.7.5

Flow rate: 1.0 ml/min

Temperature: 25°C

Detection: UV (280 nm)


Fig. 8 Separation of protein A by reversed-phase chromatography on TSKgel Octadecyl-NPR

Column : TSKgel Octadecyl-NPR (35×4.6 mm I.D.)

Sample : *Staphylococcus aureus* protein A

Elution : 10 min linear gradient of acetonitrile from 15 to 80% in 5 mM HClO₄

Flow rate: 1.5 ml/min

Temperature: 25°C

Detection: UV (220 nm)


Fig. 9 Separation of beef-soybean protein mixture by reversed-phase chromatography on TSKgel Octadecyl-NPR

Column : TSKgel Octadecyl-NPR (35×4.6 mm I.D.)

Sample : (a) beef protein (20 μg in 18 μl), (b) soybean protein (10 μg in 6 μl), (c) beef-soybean protein mixture (10 μg each)

Elution : 10 min linear gradient of acetonitrile from 5 to 80% in 100 mM HClO₄

Flow rate: 1.5 ml/min

Temperature: 25°C

Detection: UV (220 nm)

Arrows indicate characteristic peaks of soybean protein.


Fig. 10 Separation of food protein by reversed-phase chromatography on TSKgel Octadecyl-NPR

Column : TSKgel Octadecyl-NPR (35×4.6 mm I.D.)

Sample : (a) pork protein (extracted at 80°C, 5 μg in 63 μl), (b) extract of ham (5 μg in 56 μl)

Elution : 10 min linear gradient of acetonitrile from 5 to 80% in 100 mM HClO₄

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Table 1 Recovery of proteins with 400 ng loads

Column	Protein	Recovery (%)
DEAE-NPR	ovalbumin	103
	soybean trypsin inhibitor	98
	IgG ₁	89
SP-NPR	lysozyme	94
Octadecyl-NPR	ribonuclease A	95
	cytochrome c	102
	myoglobin	99


Fig. 11 Relationship between sample loads and peak area of eluate in the separation of proteins by ion-exchange chromatography on TSKgel DEAE and SP-NPR


Fig. 12 Relationship between sample loads and peak area of eluate in the separation of proteins by reversed-phase chromatography on TSKgel Octadecyl-NPR


Fig. 13 Separation of a protein mixture by reversed-phase chromatography on TSKgel Octadecyl-NPR

Column : TSKgel Octadecyl-NPR (35×4.6 mm I.D.)

Sample : 1. ribonuclease, 2. insulin, 3. cytochrome c, 4. lysozyme, 5. transferrin, 6. myoglobin (50 ng each)

Elution : 10 min linear gradient of acetonitrile from 15 to 80% in 5 mM HClO₄

Flow rate: 1.5 ml/min

Temperature: 25°C

Detection: UV (220 nm)


Fig. 14 Separation of hexokinase by ion-exchange chromatography on TSKgel DEAE-NPR

Sample loading : 500 ng
 Elution : 10 min linear gradient of NaCl from 0 to 0.5 M in 20 mM Tris-HCl buffer (pH 8.0)
 Flow rate: 1.5 ml/min
 Temperature: 25°C
 Detection: fluorescence (ex.; 280 nm, em.; 340 nm)
 Recovery of enzymatic activity was 80%.

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
Fig. 15 Separation of synthetic peptide by high-performance reversed-phase chromatography on TSKgel Octadecyl-NPR

Column : TSKgel Octadecyl-NPR (35×4.6 mm I.D.)
 Sample : dodecapeptide (EAEDLQVG-QVCG, 10 μg in 10 μl)
 Elution : 20 min linear gradient of acetonitrile from 0 to 80% in 0.1% TFA
 Flow rate: 1.5 ml/min
 Temperature: 25°C
 Detection: UV (215 nm)

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
Fig. 16 Separation of synthetic peptide by high-performance reversed-phase chromatography on TSKgel Octadecyl-NPR

Column : TSKgel Octadecyl-NPR
(35×4.6 mm I.D.)

Sample : triacontadipeptide (EAEDL-
QVGQVELGGPGAGSL-
QPLALEGSLQC, 12 μ g in
12 μ l)

Elution : 20 min linear gradient of
acetonitrile from 0 to 80% in
0.1% TFA

Flow rate: 1.5 ml/min

Temperature: 25°C

Detection: UV (215 nm)


Fig. 17 Separation of tryptic digests of BSA by reversed-phase chromatography on TSKgel Octadecyl-NPR

Column : TSKgel Octadecyl-NPR
(35×4.6 mm I.D.)

Sample : Tryptic digests of reduced
and S-carboxymethylated
BSA (10 μ g)

Eluent A : 1 mM sodium dodecyl sulfate
in 50 mM phosphate buffer,
pH 2.8

Eluent B : 60% acetonitrile in eluent A

Elution : 10 min linear gradient from
eluent A to eluent B

Flow rate: 1.5 ml/min

Temperature: 25°C

Detection: UV (210 nm)


Fig. 18 Separation of peptides by ion-exchange chromatography on TSKgel SP-NPR

Column : TSKgel SP-NPR (35×4.6
mm I.D.)

Sample : 1. γ -endorphin (0.25 μ g), 2.
bombesin (0.5 μ g), 3.
calcitonin (0.25 μ g), 4.
substance P (0.5 μ g), 5.
insulin (0.5 μ g), 6. β -endorphin
(1 μ g), 7. VIP (1 μ g), 8.
parathyroid hormone (0.5 μ g)

Elution : 10 min linear gradient of
sodium sulfate from 0 to 0.15
M in a mixture of 20 mM
acetate buffer (pH 3.5) and
acetonitrile (60/40 by volume)

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
Fig. 19 Separation of a mixture of oligonucleotides by ion-exchange chromatography on TSKgel DEAE-NPR

Column : TSKgel DEAE-NPR (35×4.6 mm I.D.)

Sample : a mixture of oligoadenylic acids with chain length of 1, 2, 3, 4, 5, 8, 10, 12, 16 and 20 nucleotides (0.1 μg each)

Elution : 20 min linear gradient of NaClO₄ from 0 to 0.4 M in 20 mM Tris-HClO₄ buffer (pH 9.0)

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
Fig. 20 Separation of oligonucleotides by ion-exchange chromatography on TSKgel DEAE-NPR
 Column : TSKgel DEAE-NPR (two 35×4.6 mm I.D. columns)
 Sample : hydrolysate of polyadenylic acid containing heptadecamer—ca. hectomer
 Elution : 60 min linear gradient of NaCl from 0.25 to 1 M in 20 mM Tris-HCl buffer (pH 9.0)
 Flow rate: 1.0 ml/min
 Temperature: 25°C
 Detection: UV (260 nm)


Fig. 21 Separation of oligonucleotides by ion-exchange chromatography on TSKgel DEAE-NPR
 Column : TSKgel DEAE-NPR (35×4.6 mm I.D.)
 Sample : 1. tetraadenylic acid (0.1 μg),
 2. tetracytidylic acid (0.2 μg),
 3. tetrathymidylic acid (0.1 μg) 4. tetraguanylic acid (0.4 μg)
 Elution : 20 min linear gradient of NaCl from 0 to 1 M in 20 mM 1,3-diaminopropane-HCl buffer (pH 10.5)
 Flow rate: 1.5 ml/min
 Temperature: 25°C
 Detection: UV (260 nm)


Fig. 22 Separation of synthetic oligonucleotides by ion-exchange chromatography on TSKgel DEAE-NPR
 Column : TSKgel DEAE-NPR (35×4.6 mm I.D.)
 Sample : synthetic oligonucleotide, 48 mer (d(AGCTAGTCA CGTACTGATGCAT CGGACTGATCGCATGCTAGTACGTCA), 6 μg)
 Elution : 20 min linear gradient of NaClO₄ from 0 to 0.4 M in 20 mM 1,3-diaminopropane-HCl buffer (pH 10.5)
 Flow rate: 1.5 ml/min
 Temperature: 25°C
 Detection: UV (260 nm)


Fig. 23 Separation of synthetic oligonucleotides by ion-exchange chromatography on TSKgel DEAE-NPR

Column : TSKgel DEAE-NPR (35×4.6 mm I.D.)

Sample : synthetic oligonucleotide, 17 mer (dGTAAACGACGGCCAGT), 1 μg)

Elution : 20 min linear gradient of NaClO₄ from 0 to 0.4 M in 20 mM Tris-HCl buffer (pH 9.0)

Flow rate: 1.5 ml/min

Temperature: 25°C

Detection: UV (260 nm)


Fig. 24 Separation of tRNA from *E. coli* by ion-exchange chromatography on TSKgel DEAE-NPR

Column : TSKgel DEAE-NPR (35×4.6 mm I.D.)

Sample : tRNA (*E. coli* MRE600, 4 μg)

Elution : 40 min linear gradient of NaCl from 0.35 M to 1 M in 20 mM Tris-HCl buffer (pH 7.5)

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
Fig. 25 Separation of pBR322-HaeIII digest by high performance ion-exchange chromatography on TSKgel DEAE-NPR

Column : TSKgel DEAE-NPR (35×4.6 mm I.D.)

Sample : pBR322-HaeIII digest (4.8 μg in 8 μl)

Elution : 0.1 min linear gradient of NaCl from 0.25 M to 0.45 M followed by 2.9 min linear gradient of NaCl from 0.45 to 0.5 M and 57 min linear gradient of NaCl from 0.5 to 1.0 M in 20 mM Tris-HCl (pH 9.0)

Flow rate: 1.5 ml/min

Temperature: 25°C

Detection: UV (260 nm)


Fig. 26 Separation of λ DNA-HindIII digest by ion-exchange chromatography on TSKgel DEAE-NPR

Column : TSKgel DEAE-NPR
(35 \times 4.6 mm I.D.)

Sample : λ DNA-HindIII digest (2 μ g in
4 μ l)

Elution : 10 min linear gradient of
NaCl from 0.5 M to 0.1 M in
20 mM Tris-HCl (pH 9)

Flow rate: 1.0 ml/min

Temperature: 25°C

Detection: UV absorbance at 260 nm


Fig. 27 Separation of rRNA from *E. coli* by ion-exchange chromatography on TSKgel DEAE-NPR

Column : TSKgel DEAE-NPR
(35 \times 4.6 mm I.D.)

Sample : 16 S, 23 S rRNA (*E. coli*
R-13, 0.1 unit)

Elution : 60 min linear gradient of
NaCl from 0.5 to 1 M in 20
mM Tris-HCl buffer (pH 7.5)

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