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Basic Evaluation of TSKgel Super Q-Toyoparl 650, a New ion-Exchanger for Preparative Separation of Proteins

Takao MATSUDA
Toshikuni KOGA
Katsuo KOMIYA
Yoshio KATO

The title new anion-exchanger, prepared by introducing quaternary ammonium groups onto TSKgel Toyoparl HW65, was evaluated for its basic characteristics and applications. This new ion-exchanger was found to possess high chemical and mechanical stability as well as high adsorption capacity for standard proteins. The recovery of proteins was almost quantitative. Accordingly, it should be particularly useful for the preparative separation of proteins.

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Table 2 Static Adsorption Capacity for BSA on Some Anion-Exchangers

Anion-Exchanger	BSA Adsorption Capacity (meq/ml-gel)
Super Q TOYOPEARL 650S	126
Super Q TOYOPEARL 650M	143
Super Q TOYOPEARL 650C	129
QAE-TOYOPEARL 550C	70
DEAE-TOYOPEARL 650M	25
Q-Sepharose Fast Flow	40
Fractogel EMD TMAE 650M	57

Adsorption Buffer: 50 mM Tris-HCl (pH 8.7)

Table 3 Dynamic Adsorption Capacity for Proteins on Some Anion-Exchangers

Anion-Exchanger	Adsorption Capacity (mg/ml-gel)		
	Thyroglobulin (MW:660,000)	BSA (MW:66,000)	Trypsin inhibitor (MW:21,500)
Super Q TOYOPEARL 650M	12	157	167
QAE-TOYOPEARL 550C	9	69	72
DEAE-TOYOPEARL 650M	12	40	47
Q-Sepharose Fast Flow	7	45	92
Fractogel EMD TMAE 650M	10	56	60

Adsorption Buffer: 50 mM Tris-HCl (pH 8.7)


Fig. 2 Chemical stability of Super Q-Toyoparl 650M in alkaline and acidic aqueous solution


Fig. 3 Chemical stability of Super Q-Toyoparl 650M in alkaline and acidic aqueous solution


Fig. 4 Relationship between flow-rate and pressure drop on Super Q-Toyoparl 650S, 650M and 650C
Column size: 15×1.6 cm I.D. (glass column)
Eluent: Deionized water

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
Fig. 5 Relationship between flow-rate and pressure drop on Super Q-Toyoparl 650S, 650M and 650C
Column size: 20×10.8 cm I.D. (stainless steel column)
Eluent: Deionized water

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
Fig. 6 Dependence of bed volume on pH and ionic strength of the buffer

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Table 4 Recovery of Proteins from Super Q-Toyopearl 650M

Protein	Recovery (%)
thyroglobulin	98
Ferritin	97
γ -globulin	101
Hemoglobin	89
Bovine serum albumin	99
Ovalbumin	100
β -lactoglobulin	100
Trypsin inhibitor	100
Myoglobin	100

Column size: 15×1.6 cm I.D. (glass column)

Proteins of 10 mg each were applied to the column equilibrated with 0.05 M Tris-HCl buffer of pH 8.7 at a flow rate of 2.0 ml/min. After elution for 1 min with 0.05 M Tris-HCl buffer of pH 8.7, the adsorbed proteins were washed out in the 0.05 M Tris-HCl buffer of pH 8.7 containing 0.5 M NaCl. The column effluent were collected for 12.5 min after sample injection and then the proteins in the effluents were determined spectrophotomerically at 280 nm.


Fig. 7 Chromatograms of a protein mixture obtained by ion-exchange chromatography on TSKgel Super Q-Toyopearl 650S, 650M and 650C

Column Size: 15×1.6 cm I.D. (glass column)

Sample: 1. carbonic anhydrase (1 mg), 2. transferrin (1.6 mg), 3. ovalbumin (2 mg), 4. trypsin inhibitor (2 mg)

Elution: 100 min linear gradient of NaCl from 0 to 0.5 M in 50 mM Tris-HCl buffer (pH 8.7)

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
Fig. 9 Separation of a protein mixture by ion-exchange chromatography on Super Q-Toyopearl 650M
 Column size: 15×1.6 cm I.D. (glass column)
 Sample: 1 carbonic anhydrase (20 mg), 2. transferrin (40 mg), 3. ovalbumin (50 mg), 4. trypsin inhibitor (50 mg)
 Elution: 100 min linear gradient of NaCl 0 to 0.5 M in Tris-HCl buffer (pH 8.7)
 Flow rate: 2.0 ml/min
 Detection: UV (280 nm)
 Temperature: 25°C


Fig. 8 Chromatogram of a mixture of ovalbumin and trypsin inhibitor obtained on Super Q-Toyopearl 650M
 Column Size: 7.5×0.75 cm I.D. (stainless steel column)
 Sample: a mixture of equal amount of ovalbumin and trypsin inhibitor
 Elution: 100 min linear gradient of NaCl from 0 to 0.5 M in 50 mM Tris-HCl buffer (pH 8.7)
 Flow rate: 1.0 ml/min
 Detection: UV (280 nm)
 Temperature: 25°C


Fig. 10 Separation of lipoxydase by ion-exchange chromatography on TSKgel Super Q-Toyopearl 650M
 Column size: 15×1.6 cm (glass column)
 Sample: a commercial crude sample of lipoxydase (30 mg)
 Elution: 100 min linear gradient of NaCl from 0 to 0.5 M in 50 mM Tris-HCl buffer (pH 8.7)
 Flow rate: 2.0 ml/min
 Detection: UV (280 nm)
 Temperature: 25°C


Fig. 11 Separation of β -lactoglobulin obtained by ion-exchange chromatography on TSKgel Super Q-Toyopearl 650S
 Column size: 15×1.6 cm I.D. (glass column)
 Sample loading: 50 mg
 Elution: 100 min linear gradient of NaCl from 0 to 0.3 M in 20 mM piperazine-HCl buffer (pH 6.0)
 Flow rate: 2.0 ml/min
 Detection: UV (280 nm)
 Temperature: 25°C

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