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Separation of DNA Fragments by High-Performance Ion-Exchange Chromatography on TSKgel DNA-NPR

Shigeru NAKATANI
Teruhiko TSUDA
Hiroyuki MORIYAMA

TSKgel DNA-NPR, a new ion-exchange column specially designed for the separation of nucleic acids, has been developed. The column packed with non-porous spherical resin-based anion-exchanger of 2.5 μm in diameter proved very useful for the separation of DNA fragments. TSKgel DNA-NPR revealed higher resolution in the separation of a wide range of DNA restriction fragments (50-23,000 base pairs) compared with a conventional non-porous anion-exchange column, TSKgel DEAE-NPR. Polymerase chain reaction (PCR) products could also be separated with high resolution. The PCR products were recovered in high yield from the column even at low level loading of sample. The detection limit of the PCR product was as low as 0.20 ng (S/N=3).

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
Fig. 1 Chromatograms of 1kb DNA ladder ($2\mu\text{g}$ in $2\mu\text{l}$).

The separations were performed on (A) a TSKgel DNA-NPR column ($4.6\text{mm I.D.} \times 7.5\text{cm L}$) and (B) a TSKgel DEAE-NPR column ($4.6\text{mm I.D.} \times 3.5\text{cm L}$) with a 20min linear gradient from 0.5 to 0.75M NaCl in 20mM Tris-HCl buffer (pH9.0) at a flow rate of $0.75\text{ml}/\text{min}$ and 25°C . DNA fragments were detected at 260nm with a $2\mu\text{l}$ flow cell. Numbers on the peaks are the chain lengths of DNA fragments in base pairs.

Table 1 Effect of void volume in the connecting tube on the resolution of DNA fragments

	Tube dimension			Resolution		
	I.D. (mm)	Length (cm)	Volume (μ l)	134bp/154bp	506bp/517bp	5090bp/6108bp
(A) Injector-Column	0.1	20	1.57	2.03	4.00	2.72
	0.1	40	3.14	2.02	3.88	2.80
	0.25	60	29.45	1.98	3.93	2.53
(B) Column-Detector	0.1	20	1.57	2.01	3.85	2.64
	0.1	40	3.14	1.95	3.53	2.59
	0.25	60	29.45	1.67	2.96	2.20

1kb DNA ladder (2μ g in 2μ l) was separated on a TSKgel DNA-NPR column under the conditions in Fig.1. The resolution was calculated from peak widths and elution volumes of each pair of DNA fragments. In (A) the connecting tube dimension between the column and the detector was 0.1mm I.D. \times 20cmL and in (B) the tube dimension between the injector and the column was 0.1mm I.D. \times 20cmL.


Fig. 2 Chromatograms of λ DNA/*Hind*III digest (1μ g in 4μ l).

The separations were performed on (A) a TSKgel DNA-NPR column (4.6mm I.D. \times 7.5cmL) and (B) a TSKgel DEAE-NPR column (4.6mm I.D. \times 3.5cmL) with a 10min linear gradient from 0.5 to 1.0M NaCl in 20mM Tris-HCl buffer (pH9.0) at a flow rate of 1.0ml/min and 25°C. DNA fragments were detected as in Fig.1. Numbers on the peaks are the chain lengths of DNA fragments in base pairs.

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Table 2 Effect of cell volume on the resolution of DNA fragments

Cell	Cell Volume (μ l)	Tube dimension between the column and the cell			Resolution		
		I.D. (mm)	Length(cm)	Volume(μ l)	134bp/154bp	506bp/517bp	5090bp/6108bp
Micro flow cell	2	0.1	20	1.57	1.93	3.51	2.30
Standard flow cell (without a heat sink)	10	0.1	20	1.57	1.90	3.24	2.15
Standard flow cell (without a heat sink)	10	0.1	40	3.14	1.84	3.21	2.21
Standard flow cell (with a heat sink)	10	0.25	100	49.09	1.38	2.47	1.84

1kb DNA ladder (2μ g in 2μ l) was separated on a TSKgel DNA-NPR column under the conditions in Fig.1. The resolution was calculated as in Table1.


Fig. 3 Dependence of the resolution on flow rate under the constant gradient time in the separation of DNA fragments.

1kb DNA ladder (2μ g in 2μ l) was separated on a TSKgel DNA-NPR column under the conditions in Fig.1 except that the flow rate was varied between 0.25 and 1.25ml/min. The resolution was calculated as in Table1.

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
Fig. 4 Dependence of the resolution on flow rate under the constant gradient volume in the separation of DNA fragments.

1kb DNA ladder (2μ g in 2μ l) was separated on a TSKgel DNA-NPR column under the conditions in Fig.1 except that the flow rate was varied between 0.25 and 1.25ml/min and the gradient time was varied between 12 and 60min to keep the gradient volume at 15ml. The resolution was calculated as in Table1.

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Table 3 Effect of response of the detector on the resolution of DNA fragments

Response (sec)	Resolution		
	134bp/154bp	506bp/517bp	5090bp/6108bp
0.05	2.08	4.16	2.80
1	1.71	3.08	2.40
3	1.20	2.03	1.82

1kb DNA ladder ($2\mu\text{g}$ in $2\mu\text{l}$) was separated on a TSKgel DNA-NPR column under the conditions in Fig.1. The resolution was calculated as in Table1.


Fig. 5 Dependence of the resolution on gradient steepness in the separation of DNA fragments.

1kb DNA ladder ($2\mu\text{g}$ in $2\mu\text{l}$) was separated on a TSKgel DNA-NPR column under the conditions in Fig.1 except that the gradient time was varied between 2.5 and 40min. The resolution was calculated as in Table1.


Fig. 6 Dependence of the resolution on temperature in the separation of DNA fragments.

1kb DNA ladder ($2\mu\text{g}$ in $2\mu\text{l}$) was separated on a TSKgel DNA-NPR column under the conditions in Fig.1 except that the temperature was varied between 15 and 45°C . The resolution was calculated as in Table1.

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
Fig. 7 Chromatograms of 1kb DNA ladder (2 μ g in 2 μ l).

The separations were performed on a TSKgel DNA-NPR column with a 15min linear gradient from 0.5 to 0.75M NaCl in 20mM Tris-HCl buffer, (A) pH7.0, (B) pH8.0 and (C) pH9.0 at a flow rate of 1.0ml/min and 25°C. DNA fragments were detected as in Fig.1. Numbers are the recoveries of DNA fragments estimated from the areas of the peaks eluted.


Fig. 8 Dependence of the resolution on sample load in the separation of DNA fragments.

1kb DNA ladder was separated on a TSKgel DNA-NPR column under the conditions in Fig.1 except that the sample load was varied between 1 and 16 μ g. The resolution was calculated as in Table1.

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
Fig. 9 Relationship between the sample load and the peak area in the separation of the DNA fragment.

PCR product of HBV DNA (from 0.38 to 3.78ng) was separated on a TSKgel DNA-NPR column with a 10min linear gradient from 0.5 to 0.75M NaCl in 20mM Tris-HCl buffer (pH9.0) at a flow rate of 0.5ml/min and 25°C. DNA fragment was detected at 260nm with a 10 μl flow cell (without a heat sink).

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
Fig. 10 Chromatogram of pBR322/HaeIII digest (2 μg in 3 μl).

The separation was performed on a TSKgel DNA-NPR column with a 0.13min linear gradient from 0.25 to 0.45M NaCl followed by a 3.87min linear gradient from 0.45 to 0.5M and a 38min linear gradient from 0.5 to 0.75M in 20mM Tris-HCl buffer (pH9.0) at a flow rate of 0.75ml/min and 35°C. DNA fragments were detected as in Fig.1. Numbers on the peaks are the estimated chain lengths of DNA fragments in base pairs.


Fig. 11 Chromatograms of PCR product of HCV DNA (290bp).

The separation was performed on a TSKgel DNA-NPR column with (A) a 20min linear gradient from 0 to 1.0M NaCl and (B) a 10min linear gradient from 0.5 to 1.0M NaCl in 20mM Tris-HCl buffer (pH9.0) at a flow rate of 0.5ml/min and 25°C. DNA fragments were detected at 260nm with a 10 μ l flow cell (without a heat sink).

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
Fig. 12 Chromatogram of PCR product of Interleukin-1 β DNA (306bp).

The separation was performed on a TSKgel DNA-NPR column (4.6mm I.D. \times 7.5cmL) with a 10min linear gradient from 0.5 to 0.75M NaCl in 20mM Tris-HCl buffer (pH9.0) at a flow rate of 0.5ml/min and 25°C. DNA fragments were detected as in Fig. 11.

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