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## The Improved Determination of Catecholamines in Plasma with Automated Catecholamine Analyzer "HLC-725CA"

Yuji HIROWATARI  
Motohiro TUBOI  
Masami TOGASHI

Approximately three percents of plasma samples for catecholamines assay in Japan cannot be determined the amount of catecholamines by the regular analysis method with HLC-725CA. The analysis method of HLC-725CA consists of two pre-treatment columns and one separation column. We considered that therapeutic drugs and their metabolites in those samples interfered with separation of catecholamines. Adequate separation of catecholamines and interference substances has been achieved by using a reversed-phase column as the separation column instead of a cation-exchange column in the regular method.

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
Fig. 2 Schematic diagram of the column-switching system of HLC-725CA

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Table 1 Reproducibility of 10 consecutive catecholamine measurements in plasma sample by the improved method

	N	Mean (ng/ml)	S.D. (ng/ml)	C.V. (%)
Epinephrine	10	0.207	0.006	2.7
Norepinephrine	10	0.777	0.022	2.9
Dopamine	10	0.167	0.003	2.2


Fig. 3 Chromatograms of plasma samples

Sample 1; A and B, sample 2; C and D, sample 3; E and F. A, C and E, and B, D, and F are analyzed patterns by the improved method, and the normal method, respectively. R.T.; retention time.


Fig. 4 Correlation of catecholamine concentrations measured by the normal method and the improved method for 39 plasma samples

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- 7) H. Nohhta, M.-K. Lee and Y. Ohkura ; Fluorescent products of the reaction for the determination of catecholamines with 1,2-diphenylethylenediamine., Anal. Chim. Acta, 267, 137-139, (1992).

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