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TSK-GEL MultiporeH Series for Gel Permeation Chromatography — Evaluation and Application —

Kuniyuki TOKUNAGA
Hiroyuki MORIYAMA

Size-exclusion chromatography (SEC) has been used widely and routinely for determining molecular mass distribution and mass average of polymers. Commonly used SEC columns are packed with a mixture of packing materials of different pore size distribution, so the chromatograms of polymers obtained on these columns often show inflectional patterns, making it difficult to determine molecular mass distribution with good reproducibility.

SEC analysis of cationic polymers using tetrahydrofuran (THF) or chloroform as the eluent is infeasible due to the interaction between the packing materials and the polymers. Now, we have developed three new types of columns packed with porous polystyrene gels of a broad pore size distribution, *i.e.*, TSKgel MultiporeH_{XL}-M, TSKgel MultiporeH_{HR}-M HT and TSKgel MultiporeH_{XL}-M C. TSK-GEL MultiporeH series show excellent linearity in calibration and no inflection is observed on the SEC chromatograms obtained. Separation of cationic polymers has been effected successfully by SEC on the column of TSKgel MultiporeH_{XL}-M C using THF as the solvent. In this paper, the characteristics as well as several SEC applications of these new types of packing materials, TSK-GEL MultiporeH series, have been presented.

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
Fig. 1 Schematic representation of calibration curves of conventional SEC packings.

- (A) Large pore size packings
- (B) Middle pore size packings
- (C) Small pore size packings




Fig. 2 Schematic representation of TSK-GEL MultiporeH packings.

- (A) Large pore size packings
- (B) Middle pore size packings
- (C) Small pore size packings

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Table 1 Characteristics of TSK-GEL MultiporeH series

Grade	Exclusion Limit (Polystyrene)	Theoretical* ¹ Plate (TP/30cm)	Particle Size (μm)	Shipping Solvent
TSKgel MultiporeH _{XL} -M	2 × 10 ⁶	16,000	5	THF* ²
TSKgel MultiporeH _{HR} -M HT	2 × 10 ⁶	16,000	5	ODCB* ³
TSKgel MultiporeH _{XL} -M C	2 × 10 ⁶	16,000	5	THF* ²

* 1 : Conditions

Column size : 7.8mmID*30cm
 Eluent : THF
 Flow rate : 1.0ml/min
 Temperature : 25℃
 Detection : UV 254nm
 Sample : Dicyclohexyl Phthalate

* 2 : Tetrahydrofuran

* 3 : o-Dichlorobenzene

Table 2 Solvent exchangeability of TSK-GEL MultiporeH series

Grade	Solvent
TSKgel MultiporeH _{XL} -M	Toluene, Xylene, Chloroform, Benzene Dichloromethane, Dichloroethane
TSKgel MultiporeH _{HR} -M HT	Chlorobenzene, Trichlorobenzene
TSKgel MultiporeH _{XL} -M C	Toluene, Xylene, Chloroform, Benzene Dichloromethane, Dichloroethane


Fig. 6 Calibration curves of TSK-GEL MultiporeH series.

Conditions

Column : (A) TSKgel MultiporeH_{XL}-M
(B) TSKgel MultiporeH_{HR}-M HT
(C) TSKgel MultiporeH_{XL}-M C

Column size : 7.8mmID*30cm

Eluent : (A)&(C) THF (B) ODCB

Flow rate : 1.0ml/min

Temperature : (A)&(C) 25°C (B) 135°C

Detection : (A)&(C) UV 254nm (B) RI

Sample : Polystyrene standards


Fig. 7 Calibration curves of TSKgel MultiporeH_{XL}-M and commercial MIXED columns.

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
Fig. 8 Calibration curves of TSKgel MultiporeH_{XL}-M and TSKgel MultiporeH_{XL}-M C with polystyrene in THF and CHCl₃.

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Table 3 Properties of TSK-GEL MultiporeH series and commercial columns

Column designation	Exclusion Limit (Polystyrene)	Mean Pore size $\Phi_m(\text{\AA})^{*1}$	Porosity (ϵ) ^{*2}	Calibration Linearity (r^2) ^{*3}	Slope (b) ^{*3}
TSKgel MultiporeH _{XL} -M	2×10^6	135	0.70	0.9993	-0.701
TSKgel MultiporeH _{HR} -M HT	2×10^6	135	0.70	0.9992	-0.753
TSKgel MultiporeH _{XL} -M C	2×10^6	135	0.70	0.9993	-0.730
MIX-BED Columns					
TSKgel GMH _{HR} -M	2.5×10^6	140	0.64	0.9991	-0.839
Commercial column (A)	4×10^6	160	0.60	0.9990	-0.932
Commercial column (B)	3×10^6	190	0.64	0.9981	-0.746

* 1 : The diameter of the probe molecules (Φ , in \AA) was calculated from the equation given by Halász and Martin²¹⁾

$$\Phi_m(\text{\AA}) = 0.62(Mw_{1/2})^{0.59}$$

Where $Mw_{1/2}$ is the gram-molecular weight of the probe

* 2 : $\epsilon = V_t - V_o / V_{col} - V_o$

Where V_t is the total permeation volume, V_o is the void volume and V_{col} is the empty column volume

* 3 : Linearity(r^2) and slope(b) of calibration were calculated from the following equations

$$r^2 = \frac{a \cdot \Sigma Y + b \Sigma XY - 1/n(\Sigma Y)^2}{\Sigma Y^2 - 1/n(\Sigma Y)^2}$$

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
Fig. 9 The effect of flow rate on theoretical plate on TSKgel MultiporeH_{XL}-M and TSKgel GMH_{HR}-M.

Conditions

Column size : 7.8mmID*30cm
 Eluent : THF
 Flow rate : 0.2~1.8ml/min
 Temperature : 25°C
 Detection : UV 254nm
 Sample : Dicyclohexyl Phthalate 20 μl(0.1%)


Fig. 11 Relationship between resolution and injection volume on TSKgel MultiporeH_{XL}-M.

Conditions

Column size : 7.8mmID*30cm
 Eluent : THF
 Flow rate : 1.0ml/min
 Temperature : 25°C
 Detection : UV 254nm
 Sample : Polystyrene standards mixture 2~100 μl
 Sample concentration : F-128(0.02%)/F-40(0.05%)
 F-10(0.05%)/F-2(0.1%)

← Inflections of the chromatograms


Fig. 10 Relationship between resolution and flow rate on TSKgel MultiporeH_{XL}-M.

Conditions

Column size : 7.8mmID*30cm
 Eluent : THF
 Flow rate : 0.2~1.8ml/min
 Temperature : 25°C
 Detection : UV 254nm
 Sample : Polystyrene standards mixture 20 μl
 Sample concentration : F-128(0.02%)/F-40(0.05%)
 F-10(0.05%)/F-2(0.1%)


Fig. 12 Chromatograms of epoxy resin on TSKgel MultiporeH_{XL}-M and conventional columns.

Conditions

Column : (A) TSKgel MultiporeH_{XL}-M
 (B) TSKgel GMH_{HR}-M
 (C) TSKgel G4000H_{XL}+G3000H_{XL}
 +G2500H_{XL}+G2000H_{XL}
 Column size : 7.8mmID*30cmx4
 Eluent : THF
 Flow rate : 1.0ml/min
 Temperature : 40°C
 Detection : RI
 Instrument : HLC-8120GPC(TOSOH)
 Sample : Epoxy resin (Epikote 1009)