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Ion Chromatography of Anions by Conductance Detection Using Sodium Naphthalenesulfonates as an Eluent

Yoshimitsu TADA
Noriko OHMURA
Tadashi OKADA

The addition of a small amount of disodium 1,5-naphthalenedisulfonate or trisodium 1,3,6-naphthalenetrisulfonate as an anti-retardant to a borate buffer eluent has been found to improve considerably the separation of inorganic anions in ion chromatography based on the conductivity detection. Ten inorganic anions including hydrophobic ones (F^- , Cl^- , NO_2^- , Br^- , NO_3^- , I^- , HPO_4^{2-} , SO_4^{2-} , SCN^- , $S_2O_3^{2-}$) were thoroughly separated within 28 minutes on a commercially available column (TSKgel IC-Anion-PW_{XL} PEEK, 4.6mm I.D. × 75mm). A good linearity has been observed for the calibration curves for six anions (F^- , Cl^- , Br^- , NO_3^- , HPO_4^{2-} , SO_4^{2-}) in the range of 1–100 ppm with relation coefficients being above 0.999. This procedure has been successfully applied to the determination of inorganic anions in tap waters and rain; chloride, nitrate, and sulfate ions were separated each other with high resolution, the relative standard deviation of reproducibility being within 1.7%.

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
Fig. 1 Effect of 1,3,6-NTS on the retention time of inorganic anions.

[Chromatographic conditions]

Column: TSKgel IC-Anion-PW_{XL} PEEK
(4.6mm I.D.×75mm)

Eluent: 1.105mM $Na_2B_4O_7$ + 25.5mM H_3BO_3

Injection volume: 100 μ l


Fig. 2 Typical chromatograms of standard inorganic anions using an alkaline buffer containing 1, 5-NDS.

Column : TSKgel IC-Anion-PW_{XL} PEEK (4.6mm I.D. × 75mm)

Eluent : (a) 39.2 μM 1, 5-NDS + 1.105mM Na₂B₄O₇ + 25.5mM H₃BO₃ (pH=7.7)

(b) 39.2 μM 1, 5-NDS + 0.1mM Na₂CO₃ + 0.4mM NaOH (pH=10.6)

(c) 46.1 μM 1, 5-NDS + 0.4mM Na₂CO₃ + 0.1mM NaOH (pH=10.2)

Injection volume : 100 μl

Peaks 1 : F⁻, 2 : Cl⁻, 3 : NO₂⁻, 4 : Br⁻, 5 : NO₃⁻, 6 : I⁻, 7 : HPO₄²⁻, 8 : SO₄²⁻,
9 : SCN⁻, 10 : S₂O₃²⁻


Fig. 3 Typical chromatograms of standard inorganic anions using a commercially available eluent (TSK eluent IC-Anion-A) and borate buffers containing 1, 3, 6-NTS.

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
Fig. 4 Calibration curves of inorganic anions using a borate buffer containing 1, 5-NDS.

Chromatographic conditions were the same as those in Fig. 2 (a).


Fig. 5 Calibration curves of inorganic anions using a borate buffer containing 1, 3, 6-NTS.

Chromatographic conditions were the same as those in Fig. 3 (b).

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
Fig. 6 Typical chromatogram of rain water.

Chromatographic conditions were the same as those in Fig. 2 (a).


Fig. 7 Typical chromatogram of tap water.

Injection volume: 50 μ l

Other chromatographic conditions were the same as those in Fig. 3 (b).

Table 1 Reproducibility for determination of tap water

Tap water (n=10)	Cl ⁻	NO ₃ ⁻	SO ₄ ²⁻
Concentration (ppm)	6.0	1.4	5.1
Peak area (mV·sec)	193.63	20.37	119.95
Standard deviation (mV·sec)	0.52	0.34	1.87
Relative standard deviation (%)	0.27	1.67	1.56

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