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Development of the Automatic Catecholamine Analyzer HLC-725CA II

Yuji HIROWATARI
Yoshimasa ITO
Masanobu KASAI
Hiroaki TAKAHASHI
Hidechika HAYASHI

The new automatic catecholamine analyzer HLC-725CA II has been developed, which can determine the amounts of three catecholamines (norepinephrine, epinephrine and dopamine) in plasma and urine samples. The instrument features a column-switching system composed of two pretreatment columns and one separation column and a high-sensitive detection unit based on a post-column reaction using a fluorogenic reagent, 1,2-diphenylethyleneamine. This system is not influenced by the presence of therapeutic drugs and their metabolites in plasma and urine samples. The detection limits of norepinephrine, epinephrine and dopamine were 14, 8 and 11 fmol/mL, respectively, under the conditions of a signal-to-noise ratio 2.

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
Fig. 3 Structure of the xenon lamp housing.

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
Fig. 4 Schematic diagram of the column-switching system of HLC-725CA II.

- : catecholamines
- ▲ : water soluble components
- : anionic components
- ★ : high hydrophobic components
- DPE : 1,2-diphenylethylenediamine
- CA : catecholamines
- MV : 6-way motor valve


Fig. 5 Chromatograms produced by HLC-725CA II.

A : standard solution (1pmol/mL)
 B : plasma for a healthy person
 C : urine for a healthy person
 R.T. : retention time

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
Fig. 6 Relationship between the peak height of the fluorescence intensity and the amount of catecholamines in standard solutions.

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
Fig. 7 Correlation of catecholamine concentrations measured by HLC-725CA and HLC-725CAII for 56 plasma samples.

NE : norepinephrine
E : epinephrine
DA : dopamine

Table 1 Reproducibility of 10 consecutive catecholamine measurements by HLC-725CAII

		N	Mean (pmol/mL)	S.D. (pmol/mL)	C.V. (%)
Norepinephrine	Standard solution	10	0.99	0.00	0.34
	Plasma	10	3.27	0.01	0.17
	Urine	10	81.65	0.31	0.38
Epinephrine	Standard solution	10	0.99	0.00	0.38
	Plasma	10	0.38	0.00	0.62
	Urine	10	15.16	0.11	0.71
Dopamine	Standard solution	10	0.99	0.01	0.66
	Plasma	10	0.16	0.00	1.20
	Urine	10	609.17	2.23	0.37


Fig. 8 Correlation of catecholamine concentrations measured by HLC-725CA and HLC-725CAII for 63 urine samples.

NE : norepinephrine
 E : epinephrine
 DA : dopamine


Fig. 9 Chromatograms of patient samples produced by HLC-725CA and HLC-725CA II.

- A: pattern of plasma for patient produced by HLC-725CA II
 B: pattern of plasma for patient produced by HLC-725CA
 C: pattern of urine for patient produced by HLC-725CA II
 D: pattern of urine for patient produced by HLC-725CA
 R.T.: retention time

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