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Development and Application of New High Performance SEC Semimicrocolumns for Water and Polar Organic Solvents

Yoshimi HASHIMOTO
Toshinao IWAEDA
Kouji SHIIGI
Hiroshi TOMIZAWA
Hiroyuki MORIYAMA

New high performance SEC semimicrocolumns for water and polar organic solvents were developed. The new columns could be used in the wide polar range solvents, namely from water to polar organic solvent, due to their excellent solvent compatibility. We investigated their fundamental properties and optimized the analytical conditions for various samples on the columns. As the results, the new concept SEC system which does not need column exchange attended with solvent exchange was established. Furthermore, the columns have higher efficiencies and smaller scale housings compared with conventional SEC columns, consequently achievement of high speed and solvent saving SEC system was suggested.

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Table 1 Specifications of TSK-GEL SuperAW series

	Exclusion Limit	Particle size	Column size
TSKgel SuperAW2500	Low	4 μ m	6.0mmID×150mm
TSKgel SuperAW3000		4 μ m	6.0mmID×150mm
TSKgel SuperAW4000	↑ ↓	6 μ m	6.0mmID×150mm
TSKgel SuperAW5000		7 μ m	6.0mmID×150mm
TSKgel SuperAW6000	High	9 μ m	6.0mmID×150mm
TSKgel SuperAWM-H	(Mixed bed column)	9 μ m	6.0mmID×150mm

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
Fig. 1 Solvent compatibility of TSK-GEL SuperAW series.

<Theoretical plates measuring conditions>

Eluent : H₂O

Flow rate : 0.6ml/min

Sample : ethylene glycol 2.5mg/ml

Injection volume : 5 μl


Fig. 2 Calibration curves of TSK-GEL SuperAW series in H₂O(1).

<Conditions>

Column size:6.0mmID×150mm

Eluent:H₂O

Flow rate:0.6ml/min

Temperature:ambient

Samples:polyethylene oxide and polyethylene glycol standards


Fig. 3 Calibration curves of TSK-GEL SuperAW series in H₂O (2).

<Conditions>

Column size:6.0mmID×150mm

Eluent:H₂O

Flow rate:0.6ml/min

Temperature:ambient

Samples:pullulan standards and oligomaltoses


Fig. 4 Calibration curves of TSK-GEL SuperAW series in MeOH.

<Conditions>

Column size:6.0mmID×150mm

Eluent:10mM LiBr in MeOH

Flow rate:0.6ml/min

Temperature:40°C

Samples:polyethylene oxide and polyethylene glycol standards


Fig. 5 Calibration curves of TSK-GEL SuperAW series in DMF.

<Conditions>

Column size:6.0mmID×150mm

Eluent:10mM LiBr in DMF

Flow rate:0.6ml/min

Temperature:40°C

Samples:polyethylene oxide and polyethylene glycol standards


Fig. 6 Calibration curves of TSK-GEL SuperAW series in DMSO.

<Conditions>

Column size:6.0mmID×150mm

Eluent:10mM NaNO₃ in DMSO

Flow rate:0.3ml/min

Temperature:40°C

Samples:pullulan standards and oligomaltoses

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
Fig. 7 Effect of additional CH₃CN concentration on calibration curve of TSKgel SuperAW2500 in H₂O.

<Conditions>

Column size:6.0mmID×150mm

Eluent:CH₃CN/H₂O

Flow rate:0.6ml/min

Temperature:ambient

Samples:polyethylene oxide and polyethylene glycol standards

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
Fig. 8 Comparison of separations of oligomers on TSKgel SuperAW2500 and TSKgel G2500PW_{XL}.

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
Fig. 13 Correlation between TSKgel SuperAWM-H and TSKgel GMPW_{XL} in the measuring of weight-average molecular weight (Mw) of polysaccharides.

<Conditions>

Column: TSKgel

SuperAWM-H (6.0mmID×150mm×2)

TSKgel GMPW_{XL} (7.8mmID×300mm×2)

Eluent: 0.2M NaNO₃

Flow rate: 0.6ml/min (TSKgel SuperAWM-H)

1.0ml/min (TSKgel GMPW_{XL})

Temperature: 40°C

Samples: sodium chondroitin sulfate (x), sodium alginate (filled diamonds), carboxymethyl cellulose (filled triangles), sodium hyaluronate (open squares) each 1 mg/ml

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
Fig. 14 Correlation between TSKgel SuperAWM-H and TSKgel α -M in the measuring of weight-average molecular weight (Mw) of DMF soluble samples.

<Conditions>

Column: TSKgel SuperAWM-H (6.0mmID \times 150mm \times 2)
TSKgel α -M (7.8mmID \times 300mm \times 2)

Eluent: 10mM LiBr in DMF

Flow rate: 0.6ml/min (TSKgel SuperAWM-H)
1.0ml/min (TSKgel α -M)

Temperature: 40 $^{\circ}$ C

Samples: ethylcellulose (○), phenol resin (resol type) (◇),
phenol resin (novolac type) (×),
vinylidene chloride-acrylonitrile copolymer (△),
polyacrylonitrile (◆), polyisopropylacrylamide (□)
each 1 mg/ml

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Temperature: 40°C

Detection: RI

Sample: ethyl hydroxyethyl cellulose 0.1mg/ml

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