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Development of a Test Reagent for Nucleic Acid Based Infectious Disease. Detection of *Mycobacterium Tuberculosis*

Noriyoshi MASUDA
 Shigeo TSUCHIYA
 Yuichi ISAWA
 Takumi TOKUNAGA
 Masahide KONDO
 Ryuichi HORIE
 Kiyoshi YASUKAWA
 Takahiko ISHIGURO

We have recently developed the TRC method that is based on a homogeneous detection system involving an INAF (intercalation-activating fluorescence) DNA probe and an isothermal RNA amplification through the TRC (transcription-reverse transcription concerted) reaction. In this paper, we report an effective test reagent for nucleic acid based infectious disease, i.e., a specific reagent to detect *Mycobacterium tuberculosis*. The target RNA of the reagent is mRNA of *Pab* DNA encoding the protein antigen b specific to *Mycobacterium tuberculosis*. Time course measurement of the fluorescence intensity at 520 nm of the reaction mixture was carried out using our trial product for real-time monitoring. Since the intensity of fluorescence was proportional to the initial amount of *Pab* standard RNA (10^2 - 10^8 copies/test), the amount of *Pab* mRNA was readily determinable from the detection time. Fluorescence enhancement was specific to the nucleic acids extracted from the sputum of the patients infected with pulmonary tuberculosis, but not so to those from the ones infected with *Mycobacterium avium* complex (MAC). The BCG cells were cultured in the presence of anti-tuberculosis drugs and after 24 h the decrease in amount of *Pab* mRNA was confirmed by the TRC method using this reagent. This results demonstrate that the TRC method is useful for detecting the viable cells and can be applied to the monitoring of the patients infected with pulmonary tuberculosis and also to the drug susceptibility test of *Mycobacterium tuberculosis*.

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
Fig. 1 TRC protocol. Fluorescence intensity of the reaction mixture is measured using an instrument named TRC monitor. The reagent and enzyme mixtures were added to a sample in this order.


Fig. 2 Profile of fluorescence enhancement of starting copies of *Pab* standard RNA. A broad line shows the cut-off value corresponding to a fluorescence intensity ratio of 1.2.

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
Fig. 3 Calibration curve. Detection time, which corresponds to a time reaching the cut-off value, is plotted against the logarithms of starting copies of *Pab* RNA (10^2 - 10^8 copies). Error bars mean \pm SD.

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
Fig. 4 Detection of *Pab* mRNA in nucleic acids extracted from the sputum of the patients infected with pulmonary tuberculosis, pulmonary *Mycobacterium avium* complex (MAC) and others.

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
Fig. 5 *Pab* mRNA from BCG cultured in the presence of drugs. The BCG cells were cultured in the presence or absence of an anti-tuberculosis drug. The amount of *Pab* mRNA in 1 mL of the culture liquid is determined by the TRC method and shown as the relative values against those obtained in the absence of drugs.


Fig. 6 Electrophoresis of PCR products. *Pab* DNA was amplified by PCR from the nucleic acids extracted from BCG cells, cultured in the presence or absence of kanamycin (20 μ g/mL) for a given hour, and then subjected to 1.5% agarose gel treatment. The arrow indicates the amplified *Pab* DNA.

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